




# TruSight<sup>®</sup> Cancer Workflow on the MiniSeq<sup>™</sup> System



This checklist supports the following workflow choices.

Sequencing Instrument:	MiniSeq System
Setup Option	Local Run Manager
Library Preparation Kit:	TruSight Rapid Capture
Indexing:	Dual Indexing
Reagent Kits:	MiniSeq Kit
Analysis Workflow:	DNA Enrichment
Analysis Software:	Local Run Manager

MiniSeq sequencing using TruSight Rapid Capture libraries\_

## Set Run Parameters

- 1 Log in to Local Run Manager.
- 2 Click **Create Run**, and select **DNA Enrichment**.
- 3 Enter a run name that identifies the run.
- 4 [Optional] Enter a run description.
- 5 From the Library Kit drop-down list, select TruSight Enrichment Panels.
- 6 Specify the number of index reads.
- 7 Specify a read type: **Single Read** or **Paired End**.
- 8 Enter the number of cycles for the run.
- 9 Select an alignment method.
- 10 Select a variant calling method.
- 11 Specify the **Manifest Padding** threshold.
- 12 Enable or disable the Flag PCR Duplicates and Indel Realignment settings.
- 13 Click **Show advanced module settings** and enable or disable Picard HS metrics.
- 14 Click **Import Manifests**.
- 15 Navigate to the manifest file.
- 16 Enter a unique sample ID.
- 17 [Optional] Enter a sample description.
- 18 Select an Index 1 adapter.
- 19 Select an Index 2 adapter.
- 20 Select a manifest file.
- 21 Click **Save Run**.

## Tagment Genomic DNA

- 1 Quantify gDNA using a fluorometric method.
- 2 Dilute gDNA in Tris-HCl 10 mM, pH 8.5 to a final volume of 10 µl at 5 ng/µl.
- 3 Add the following to a new plate.
  - ▶ Normalized gDNA (10 µl)
  - ▶ TD (25 µl)
  - ▶ TDE1 (15 µl)
- 4 Shake at 1800 rpm for 1 minute.
- 5 Centrifuge at 280 × g for 1 minute.
- 6 Place on the 58°C microheating system with the lid closed for 10 minutes.
- 7 Add 15 µl ST.
- 8 Shake at 1800 rpm for 1 minute.
- 9 Centrifuge at 280 × g for 1 minute.
- 10 Incubate at room temperature for 4 minutes.

## Clean Up Tagmented DNA

- 1 Add 65 µl SPB.
- 2 Shake at 1800 rpm for 1 minute.
- 3 Incubate at room temperature for 8 minutes.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on a magnetic stand until liquid is clear.
- 6 Remove and discard all supernatant.
- 7 Wash 2 times with 200 µl 80% EtOH.
- 8 Use a 20 µl pipette to remove residual EtOH.
- 9 Air-dry on the magnetic stand for 10 minutes.
- 10 Remove from the magnetic stand.
- 11 Add 22.5 µl RSB.
- 12 Shake at 1800 rpm for 1 minute.
- 13 Incubate at room temperature for 2 minutes.
- 14 Centrifuge at 280 × g for 1 minute.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Transfer 20 µl supernatant.

## Amplify Tagmented DNA

- 1 Arrange Index 1 (i7) adapters in columns 1–12.
- 2 Arrange Index 2 (i5) adapters in rows A–H.
- 3 Place the plate on the TruSeq Index Plate Fixture.
- 4 Add 5 µl of each Index 1 adapter down each column.
- 5 Add 5 µl of each Index 2 adapter across each row.
- 6 Add 20 µl NLM.
- 7 Shake at 1200 rpm for 1 minute.
- 8 Centrifuge at 280 × g for 1 minute.
- 9 Place on the thermal cycler and run the NLM AMP program.

**SAFE STOPPING POINT**

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

## Clean Up Amplified DNA

- 1 Centrifuge at 280 × g for 1 minute.
- 2 Transfer 50 µl supernatant.
- 3 Add 90 µl SPB.
- 4 Shake at 1800 rpm for 1 minute.
- 5 Incubate at room temperature for 10 minutes.
- 6 Centrifuge at 280 × g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200 µl 80% EtOH.
- 10 Use a 20 µl pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 10 minutes.
- 12 Add 27.5 µl RSB.
- 13 Shake at 1800 rpm for 1 minute.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at 280 × g for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 25 µl supernatant.
- 18 Quantify the library using a fluorometric method.

**SAFE STOPPING POINT**

If you are stopping, seal the plate and store at -25°C to -15°C for up to 14 days.

## Hybridize Probes

- 1 Combine 500 ng of each DNA library. Make sure that each library has a unique index.
  - ▶ For total volume > 40  $\mu$ l, concentrate the pooled sample to 40  $\mu$ l.
  - ▶ For total volume < 40  $\mu$ l, increase the volume to 40  $\mu$ l with RSB.
- 2 Use 500 ng of each DNA library quantified by QuantiFluor.
- 3 Add the following to a new plate.
  - ▶ Library pool (40  $\mu$ l)
  - ▶ EHB (50  $\mu$ l)
  - ▶ CSO (10  $\mu$ l)
- 4 Shake at 1200 rpm for 1 minute.
- 5 Centrifuge at 280  $\times$  g for 1 minute.
- 6 Place on the thermal cycler and run the NRC HYB program.
- 7 Keep at the 58°C holding temperature for at least 90 minutes and up to 24 hours.

## Capture Hybridized Probes

- 1 Centrifuge at 280  $\times$  g for 1 minute.
- 2 Transfer all volumes.
- 3 Add 250  $\mu$ l SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at 280  $\times$  g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash 2 times with 200  $\mu$ l EWS.
- 11 Mix 28.5  $\mu$ l EE1 and 1.5  $\mu$ l 2 N NaOH, and then vortex.
- 12 Add 23.5  $\mu$ l elution premix.
- 13 Shake at 1800 rpm for 2 minutes.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at 280  $\times$  g for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 21  $\mu$ l supernatant.
- 18 Add 4  $\mu$ l ET2.
- 19 Shake at 1200 rpm for 1 minute.
- 20 Centrifuge at 280  $\times$  g for 1 minute.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

## Perform Second Hybridization

- 1 Add the following.
  - ▶ RSB (15  $\mu$ l)
  - ▶ EHB (50  $\mu$ l)
  - ▶ CSO (10  $\mu$ l)
- 2 Shake at 1200 rpm for 1 minute.
- 3 Centrifuge at 280  $\times$  g for 1 minute.
- 4 Place on the thermal cycler and run the NRC HYB program.
- 5 Keep at the 58°C holding temperature for at least 14.5 hours and up to 24 hours.

## Perform Second Capture

- 1 Centrifuge at  $280 \times g$  for 1 minute.
- 2 Transfer supernatant.
- 3 Add 250  $\mu$ l SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash 2 times with 200  $\mu$ l EWS.
- 11 Mix 28.5  $\mu$ l EE1 and 1.5  $\mu$ l 2 N NaOH, and then vortex.
- 12 Add 23.5  $\mu$ l elution premix.
- 13 Shake at 1800 rpm for 2 minutes.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at  $280 \times g$  for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 21  $\mu$ l supernatant.
- 18 Add 4  $\mu$ l ET2.
- 19 Shake at 1800 rpm for 1 minute.
- 20 Centrifuge at  $280 \times g$  for 1 minute.

## Clean Up Captured Library

- 1 Add 45  $\mu$ l SPB.
- 2 Shake at 1800 rpm for 1 minute.
- 3 Incubate at room temperature for 10 minutes.
- 4 Centrifuge at  $280 \times g$  for 1 minute.
- 5 Place on a magnetic stand until liquid is clear.
- 6 Remove and discard all supernatant.
- 7 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 8 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 9 Air-dry for 10 minutes.
- 10 Add 27.5  $\mu$ l RSB.
- 11 Shake at 1800 rpm for 1 minute.
- 12 Incubate at room temperature for 2 minutes.
- 13 Centrifuge at  $280 \times g$  for 1 minute.
- 14 Place on a magnetic stand until liquid is clear.
- 15 Transfer 25  $\mu$ l supernatant.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

## Amplify Enriched Library

- 1 Add 5  $\mu$ l PPC.
- 2 Add 20  $\mu$ l NEM.
- 3 Shake at 1200 rpm for 1 minute.
- 4 Centrifuge at  $280 \times g$  for 1 minute.
- 5 Place on the thermal cycler and run the NEM AMP12 program.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $2^{\circ}\text{C}$  to  $8^{\circ}\text{C}$  for up to 2 days.

## Clean Up Amplified Enriched Library

- 1 Centrifuge at  $280 \times g$  for 1 minute.
- 2 Transfer 50  $\mu$ l.
- 3 Add 90  $\mu$ l SPB.
- 4 Shake at 1800 rpm for 1 minute.
- 5 Incubate at room temperature for 10 minutes.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 10 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 10 minutes.
- 12 Add 32.5  $\mu$ l RSB.
- 13 Shake at 1800 rpm for 1 minute.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at  $280 \times g$  for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 30  $\mu$ l supernatant.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

## Check Enriched Libraries

- 1 Quantify using a fluorometric method.
- 2 If the concentration is higher than the quantitative range for the High Sensitivity DNA chip, dilute the library 1:10 with RSB.
- 3 Run 1  $\mu$ l using a High Sensitivity DNA chip.

## Prepare Consumables

- 1 Remove the reagent cartridge from  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  storage.
- 2 Thaw reagents in a room temperature water bath for 90 minutes.
- 3 Invert the cartridge 5 times to mix reagents.
- 4 Gently tap on the bench to reduce air bubbles.
- 5 Remove a new flow cell package from  $2^{\circ}\text{C}$  to  $8^{\circ}\text{C}$  storage.
- 6 Set the unopened flow cell package aside at room temperature for 30 minutes.
- 7 Remove the flow cell from the foil package and flow cell container.
- 8 Clean the glass surface of the flow cell with a lint-free alcohol wipe.
- 9 Dry with a lint-free lens cleaning tissue.



## Denature, Dilute, and Load Libraries

- 1 Dilute 100 µl 1 N NaOH to 1 ml 0.1 N NaOH.
- 2 Invert the tube several times to mix.
- 3 Thaw the Hybridization Buffer at room temperature.
- 4 Vortex briefly before use.
- 5 Thaw the RSB at room temperature.
- 6 Transfer 25 µl of the 4 nM library pool to a new microcentrifuge tube.
- 7 Add 75 µl RSB to dilute to 1 nM.
- 8 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 9 Combine 5 µl library with 5 µl 0.1 N NaOH.
- 10 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 11 Incubate at room temperature for 5 minutes.
- 12 Add 5 µl 200 mM Tris-HCl, pH 7.0.
- 13 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 14 Add 985 µl of prechilled Hybridization Buffer.
- 15 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 16 Transfer 180 µl library to a new microcentrifuge tube.
- 17 Add 320 µl prechilled Hybridization Buffer.
- 18 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 19 [Optional] Denature and dilute a PhiX control to 1.8 pM and a 1% spike-in to the final library.
- 20 Clean the foil seal covering reservoir #16 using a low-lint tissue.
- 21 Pierce the seal with a clean 1 ml pipette tip.
- 22 Add 500 µl prepared libraries into reservoir #16.

## Perform a Sequencing Run

- 1 From the Home screen, select **Sequence**.
- 2 Enter your user name and password.
- 3 Select **Next**.
- 4 Select a run name from the list of available runs.
- 5 Select **Next**.
- 6 Open the flow cell compartment door.
- 7 Press the release button to the right of the flow cell latch.
- 8 Place the flow cell on the flow cell stage over the alignment pins.
- 9 Close the flow cell latch to secure the flow cell.
- 10 Close the flow cell compartment door.
- 11 Open the reagent compartment door.
- 12 Slide the reagent cartridge into the reagent compartment until the cartridge stops.
- 13 Remove the spent reagents bottle from the compartment.
- 14 Discard the contents and slide the empty spent reagents bottle into the compartment.
- 15 Close the compartment door and select **Next**.
- 16 Confirm run parameters.
- 17 Select **Next**.
- 18 When the automated check is complete, select **Start**.
- 19 Monitor run progress, intensities, and quality scores as metrics appear on the screen.

## View Analysis Results

- 1 From the Local Run Manager dashboard, click the run name.
- 2 From the Run Overview tab, review the sequencing run metrics.
- 3 [Optional] Click the **Copy to Clipboard**  icon for access to the output run folder.
- 4 Click the Sequencing Information tab to review run parameters and consumables information.
- 5 Click the Samples and Results tab to view the analysis report.
- 6 [Optional] Click the **Copy to Clipboard**  icon for access to the Analysis folder.