

# Nextera® Library Prep for the MiSeq® System

Sequencing's fastest library prep delivers quality *de novo* assembly of small genomes.

## Introduction

Next-generation sequencing is emerging as an important tool for both resequencing and *de novo* assembly of small genomes. An integral part of the MiSeq system's streamlined sequencing workflow is Nextera Library Preparation. With Nextera, sequencer-ready libraries can be prepared in 90 minutes with only 15 minutes of hands-on time (Figure 1), enabling researchers to go from 50 ng of sample DNA to analyzed data in as little as 8 hours with the MiSeq system. This application note shows how the unique transposome-mediated fragmentation and adapter ligation enables rapid and efficient library prep for rapid short read sequencing on the MiSeq system, as well as how Nextera library prep supports the highest quality sequence for *de novo* assembly using longer, paired-end reads.

## Results

Genomic DNA isolated from the well-characterized *Escherichia coli* strain MG1655 was used to prepare a sequencing library in 90 minutes using Illumina's Nextera DNA Library Preparation Kit. For sequencing on the MiSeq instrument, prepared samples were placed in the reagent cartridge and loaded on the instrument along with the flow cell. All subsequent steps were performed on the instrument, including cluster generation, single or paired-end sequencing, and primary data analysis.

For resequencing, automated data analysis was performed directly on the MiSeq integrated computer, requiring no specialized servers or computing facilities. Primary run metrics of the Nextera-prepped 1 × 35 bp run completed in 8 hours are shown in Table 1. Nextera library prep results in sequence data with high cluster densities, excellent quality scores, and high genome coverage.


Table 1: Primary Resequencing Metrics 1 × 35 bp

Metric	MiSeq System
Total Time*	8 hours
Passing Filter Reads (million)	5.82
% Bases ≥ Q30 (99.9% accuracy)	94.5
Mean Coverage Depth	41×
% Genome Coverage	97.89
Number of Gaps	345
Average Gap Size	296
Total Gaps (Kb)	102
Total SNPs	32

\*Total time includes library prep from genomic DNA with Nextera, 1 × 35 bp single-read sequencing run, and automated, on-instrument data analysis.

For *de novo* sequencing, the *E. coli* library prepared using Nextera library prep reagents was sequenced using a 2 × 150 read length on MiSeq. A previously-described *E. coli* strain MG1655 prepared using TruSeq® library preparation reagents was sequenced on MiSeq at 2 × 150 bp for comparison<sup>1</sup>. Both assemblies were performed using Velvet v1.1.05<sup>2</sup>. The *de novo* assembly results for samples prepared with both Nextera and TruSeq reagents are shown in Table 2. Both library prep methods deliver excellent *de novo* assembly results.

Figure 1: Nextera Delivers the Fastest Sequencer-Ready Library Prep



Nextera library preparation requires only 90 minutes with 15 minutes of hands-on time. Unlike other personal sequencing systems<sup>3</sup>, no auxiliary equipment other than a standard thermocycler is required.

