





# More informed PGD.

Two methods of PGD are currently in use: STR analysis and karyomapping (Table 1).

## STR analysis

Current PGD methods rely on examining STRs adjacent to specific disease loci to identify a gene mutation for a particular disorder. STRs are repeating sequences of 2–6 base pairs of DNA; a common STR is the sequence CA. STRs often differ in repeat number between alleles and between individuals. When used for PGD, each STR test for a genetic disorder must be developed individually, resulting in a time-consuming process that can only be accessed at a few specialist laboratories. Each test examines only a few STR markers and can produce inaccurate results due to recombination events.

To determine the source of parental chromosomes in the embryo, a laboratory performing STR analysis searches for microsatellite regions surrounding the gene of interest that contains specific STRs. If microsatellites can be identified, it is possible to use PCR to calculate the size of each of microsatellite region. If the sizes of these regions differ between each parent, it is then possible to identify the origin, maternal or paternal, of the chromosome region the child has inherited.

STR analysis has a number of inherent challenges:

1. Microsatellite regions must be on both sides of the gene of interest to make sure that a recombination event has not occurred within that region.
2. Microsatellite regions may not be informative if, for example, the parents happen to have the same number of STRs in the microsatellite region of interest.
3. Microsatellite regions can be a significant distance from the gene, reducing confidence that a recombination event has not 'swapped' the affected gene into that region without affecting the flanking microsatellite regions.
4. PCR probes must be designed for each different condition, may need to be specific to the family being studied, and may not always work.
5. Interpreting STRs is a skilled, time-consuming, and sometimes subjective, process.






# About karyomapping.

Karyomapping uses SNP genotyping data from the parents and reference to create a comprehensive map of the parental origin of chromosome segments inherited by the embryo. These chromosome segments are called haploblocks (Figure 2). By establishing the inheritance of haploblocks surrounding the region of interest, it is possible to infer the disease status of each embryo—affected, carrier, or unaffected (Figure 3)—by comparing it to the disease status of a reference.

Figure 3: Inherited Haploblocks



A pair of chromosomes in an embryo with their predicted haploblocks. The green and yellow chromosome has come from the mother, the red and blue chromosome has come from the father.

Figure 4: Identifying the Inheritance Status of an Embryo




Illustration of the inheritance possibilities for four embryos from a karyomapping case with a recessive disorder. By identifying which haploblocks have been inherited, the status of the embryo (affected, carrier, or unaffected) can be established.

# Glossary

## 1000 Genomes Project (1kGP)

An international collaboration to produce an extensive public catalog of human genetics variation, including SNPs; learn more at [www.1000genomes.org](http://www.1000genomes.org)

## Haploblock

Intact chromosome segment between recombination sites that is inherited from one of the parental chromosomes

## In vitro Fertilization (IVF)

A complex procedure that involves removal of eggs, fertilization of the eggs in a laboratory dish, and transfer of an embryo into the uterus

## Karyomapping

A comprehensive method for genome-wide linkage-based analysis of single-gene defects that can be used on a single cell or small number of cells from an embryo in preimplantation genetic diagnosis (PGD)

## Microsatellite

Simple DNA sequence repeats of 2 to 6 bases that can differ between individuals and chromosomes

## Preimplantation Genetic Diagnosis (PGD)

The genetic testing of embryos before implantation in the uterus

## Recombination

The process by which two chromosomes exchange genetic information during meiosis

## Reference

The sibling, or other closely related relative, used to identify the affected haplotype

## Short Tandem Repeat (STR) Analysis

A method for comparing specific DNA loci from two or more samples

## Single Nucleotide Polymorphism (SNP)

A DNA sequence variation that occurs at a single nucleotide

## Tag SNP

A SNP that can be used to represent a genomic region with high linkage disequilibrium

# Learn more.

Karyomapping provides a faster, easier method for PGD before embryo implantation in IVF. To learn more or to find a specialist PGD laboratory, visit [www.illumina.com/PGD](http://www.illumina.com/PGD).

## Reference

1. Handside AH, Harton GL, Mariani B, Thornhill AR, Affara N, et al. (2010) Karyomapping: A universal method for genome wide analysis of genetic disease based on mapping crossovers between parental haplotypes. *J Med Genet.* 47: 651–658.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

© 2014 Illumina, Inc. All rights reserved.  
Illumina, Genetic Energy, Infinium, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.  
Pub No. 1570-2014-001 Current as of 24 March 2014



AGAAATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTT  
TCAAAGTACCGTAACGAAACGTATCAATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTT  
CGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTT  
AACCTACCATTAAGAGTACCGTGCAACGAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTT  
AGAAATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTT  
GATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTT  
CGTATCAATTGAGACTAAATATTAACGTACCATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTTCTGTAAACCTT