

Fragment DNA

- 1 Quantify gDNA.
- 2 Normalize gDNA with RSB in a final volume of 15 μ l in a new plate:
 - ▶ For a 350 bp insert size—1.67 ng/ μ l of sample
 - ▶ For a 550 bp insert size—5 ng/ μ l of sample
- 3 Transfer 15 μ l DNA to a microTUBE-15.
- 4 Centrifuge at 3000 \times g for 1 minute.
- 5 Fragment the DNA on a Covaris.

Table 1 Covaris S220 or E220 Settings

Setting	350 bp Insert	550 bp Insert
Duty factor		20%
Peak Incident Power		18 W
Cycles per burst		50
Duration	45 seconds	22 seconds
Temperature		20°C
Water Level—S220		15
Water Level—E220		10

Table 2 Covaris M220 Settings

Setting	350 bp Insert	550 bp Insert
Duty factor		20%
Peak Incident Power		30 W
Cycles per burst		50
Duration	42 seconds	23 seconds
Temperature		20°C

- 6 Centrifuge at 600 \times g for 5 seconds.
- 7 Transfer 15 μ l DNA to a new plate.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Prepare Samples for Loading

- 1 Add 1 of the following to a new 1.5 ml microcentrifuge tube:
 - ▶ For a 350 bp insert size—900 μ l SC350
 - ▶ For a 550 bp insert size—700 μ l SC550
- 2 Vortex DMB.
- 3 Add 100 μ l DMB. Vortex for 5 seconds.
- 4 Pour the DMB and SC mixture into a new reservoir.
- 5 Add 35 μ l DMB and SC mixture to the sample plate. Pipette to mix.
- 6 Shake or vortex at 1400 rpm for 12 minutes.

Set Up Run and Load Library Card

- 1 Vortex the reagent plate for 3 seconds.
- 2 Centrifuge at 600 \times g for 5 seconds.
- 3 Select **Prepare Libraries**.
- 4 Select a protocol or run. Select **Next**.
- 5 Configure the run. Select **Next**.
- 6 Review the run and sample information. Select **Next**.
- 7 Enter tracking information. Select **Next**.
- 8 Place the library card on the stage.
- 9 Close the door. Select **Verify Library Card**.
- 10 Place the guide on the library card.
- 11 Load the oil.
- 12 Transfer 45 μ l of samples 1–8 and 9–16.
- 13 Add 45 μ l RSB to empty sample wells.
- 14 Transfer 125 μ l of reagents i–iv and v–vii.
- 15 Vortex DMB.
- 16 Add 60 μ l DMB to reagent well viii.
- 17 Transfer 15 μ l of reagents 1–4 and 5–8.
- 18 Transfer 5 μ l of reagents a–d and e–h.
- 19 Transfer 3 μ l of adapters A–H and I–P.
- 20 Remove the guide.
- 21 Close the door. Select **Start Run**.
- 22 When the run is complete, select **Next**.

Unload Libraries

- 1 Add 10 µl RSB to a new plate.
- 2 Place the guide on the library card.
- 3 Transfer 20 µl from 1L–8L and 9L–16L to the plate. Pipette to mix.
- 4 Centrifuge briefly.
- 5 Transfer from plate wells 1–8 and 9–16 to the library separation tubes.
- 6 Let stand for 10 seconds.
- 7 Transfer from library separation tubes 1–8 and 9–16 to a new plate.
- 8 Remove the library card and guide.
- 9 Select **Home**.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 2 months.

Pool Libraries

- 1 Determine the number of samples to combine.
- 2 Transfer 5 µl to a single well of a new plate. Pipette to mix.
- 3 Proceed to cluster generation.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 2 months.

Acronyms

Acronym	Definition
DMB	Digital Microfluidics Beads
RSB	Resuspension Buffer
SC350	Sample Concentration Solution 350
SC550	Sample Concentration Solution 550