


TruSight™ Oncology Comprehensive (EU)

Learn more about this CE-marked *in vitro* diagnostic, next-generation sequencing solution for comprehensive genomic profiling

As a global leader in next-generation sequencing (NGS) and microarray-based solutions, Illumina is dedicated to improving human health by unlocking the power of the genome. Illumina continues to innovate by offering TruSight Oncology Comprehensive (EU), a CE-marked *in vitro* diagnostic (IVD), pan-cancer comprehensive genomic profiling (CGP) panel. TruSight Oncology Comprehensive (EU) can generate a broad molecular profile of solid tumor patient samples, including formalin-fixed, paraffin-embedded (FFPE) tissue, maximizing a lab's ability to find actionable alterations that can help inform therapy decisions according to clinical guidelines.

With TruSight Oncology Comprehensive (EU), your laboratory can:



Generate a CGP report for a patient sample

- Detect DNA plus RNA variants and biomarkers signatures for multiple solid tumor types, generate a CGP report for a patient's tumor, and increase confidence in treatment decisions



Enable targeted therapies and clinical trials

- Leverage content that includes key biomarkers associated with drug labels, [ESMO](#) guidelines, and clinical trials for multiple solid tumor types
- Deliver results that inform therapy decisions according to clinical guidelines



Perform IVD testing in-house

- Implement a streamlined workflow, going from sample to report in 4–5 days
- Enable precision oncology testing while keeping data securely within your institution and eliminating the need for send-out services

About TruSight Oncology Comprehensive (EU)

- TruSight Oncology Comprehensive (EU) can be implemented easily in-house, features a streamlined workflow that proceeds from sample to report in 4–5 days, and requires as few as five FFPE slides
- The IVD test reliably detects all DNA and RNA variant categories, including single nucleotide variants (SNVs), insertions/deletions (indels), amplifications, fusions, and splice variants, and it enables analysis of the genomic signatures microsatellite instability (MSI) and tumor mutational burden (TMB)
- Illumina is advancing a companion diagnostic (CDx) roadmap with multiple indications linked to breakthrough therapies to improve patient outcomes

Frequently asked questions

Who is TruSight Oncology Comprehensive (EU) for?

TruSight Oncology Comprehensive (EU) is designed for molecular pathologists and lab directors who perform solid-tumor testing and provide results that guide oncologists in treatment decisions.

Types of institutions:

- Academic medical centers
- Large- to medium-sized hospitals
- Independent commercial laboratories

TruSight Oncology Comprehensive (EU) will be available in the following countries:

- | | | | |
|------------------|-----------------|---------------|------------------------|
| • Austria | • Germany | • Macedonia | • Slovenia |
| • Belgium | • Greece | • Malta | • South Africa |
| • Bulgaria | • Iceland | • Netherlands | • Spain |
| • Cyprus | • Ireland | • Norway | • Sweden |
| • Czech Republic | • Israel | • Poland | • Switzerland |
| • Denmark | • Italy | • Portugal | • United Arab Emirates |
| • Estonia | • Latvia | • Qatar | • United Kingdom |
| • Finland | • Liechtenstein | • Romania | |
| • France | • Lithuania | • Serbia | |
| | • Luxembourg | • Slovakia | |

What are the IVD claims?

TruSight Oncology Comprehensive (EU) is an IVD test that uses targeted NGS to detect variants in 517 genes using nucleic acids extracted from FFPE tumor tissue samples from patients with solid malignant neoplasms using the Illumina NextSeq™ 550Dx Instrument. The test can be used to detect SNVs, multinucleotide variants, indels, and gene amplifications from DNA, and gene fusions and splice variants from RNA. The test also reports MSI status and a TMB score.

The test is intended as a CDx to identify patients with cancer for treatment with the targeted therapy listed in [Table 1](#), in accordance with the approved therapeutic product labeling. In addition, the test is intended to provide tumor profiling information for use by qualified health care professionals in accordance with professional guidelines and is not conclusive or prescriptive for labeled use of any specific therapeutic product.

Table 1: CDx indication

Tumor type	Biomarkers	Targeted therapy
Solid tumors	<i>NTRK1, NTRK2, NTRK3</i> gene fusions	VITRAKVI® (larotrectinib)

Read the [TruSight Oncology Comprehensive Package Insert](#) to learn more.

What gene/biomarker content is included?

TruSight Oncology Comprehensive (EU) includes key biomarkers in clinical guidelines, drug labels, and clinical trials, across multiple solid tumor types and histologies ([Table 2](#)). Content includes small DNA variants ([Table 3](#)) as well as fusions (RNA), splice variants (RNA), amplifications (DNA), and complex genomic signatures ([Table 4](#)). These lists are based on content validated by Illumina.

Table 2: Subset of genomic tumor profiling biomarkers for multiple cancer types









Tumor type		Select genes with biomarkers of clinical significance ^a
	Pan-cancer	<i>BRAF, FGFR1, FGFR2, FGFR3, MSI, NTRK1, NTRK2, NTRK3, RET, TMB</i>
	Breast cancer	<i>AKT1, BRCA1, BRCA2, ERBB2, ESR1, PALB2, PIK3CA, PTEN</i>
	Colorectal cancer	<i>BRAF, ERBB2, KRAS, NRAS, POLE, MSI</i>
	Melanoma	<i>BRAF, KIT, NRAS</i>
	Non-small cell lung cancer	<i>ALK, BRAF, EGFR, ERBB2, KRAS, MET, NRG1, RET, ROS1</i>
	Ovarian cancer	<i>BRCA1, BRCA2</i>
	Pancreatic cancer	<i>BRCA1, BRCA2, KRAS, PALB2, NRG1</i>
	Prostate cancer	<i>ATM, BRCA1, BRCA2, PALB2, PTEN</i>
<p>a. Genes with biomarkers of clinical significance linked to major oncology guidelines. MSI, microsatellite instability; TMB, tumor mutational burden.</p>		

Table 3: TruSight Oncology Comprehensive (EU) panel content

Small variants: 517 genes (from DNA)											
ABL1	BCORL1	CIC	EPHB1	FGF23	HIST1H3B	JUN	MLH1	PAK5	PRKAR1A	RYBP	TAF1
ABL2	BCR	COP1	ERBB2	FGFR1	HIST1H3C	KAT6A	MLLT3	PALB2	PRKCI	SDHA	TBX3
ABRAXAS1	BIRC3	CREBBP	ERBB3	FGFR2	HIST1H3D	KDM5A	MPL	PARP1	PRKDC	SDHAF2	TCF3
ACVR1	BLM	CRKL	ERBB4	FGFR3	HIST1H3E	KDM5C	MRE11	PAX3	PRKN	SDHB	TCF7L2
ACVR1B	BMPRI1A	CRLF2	ERCC1	FGFR4	HIST1H3F	KDM6A	MSH2	PAX5	PRSS8	SDHC	TERC
ADGRA2	BRAF	CSF1R	ERCC2	FH	HIST1H3G	KDR	MSH3	PAX7	PTCH1	SDHD	TERT
AKT1	BRCA1	CSF3R	ERCC3	FLCN	HIST1H3H	KEAP1	MSH6	PAX8	PTEN	SETBP1	TET1
AKT2	BRCA2	CSNK1A1	ERCC4	FLI1	HIST1H3I	KEL	MST1	PBRM1	PTPN11	SETD2	TET2
AKT3	BRD4	CTCF	ERCC5	FLT1	HIST1H3J	KIF5B	MST1R	PDCD1	PTPRD	SH2D1A	TFE3
ALK	BRIP1	CTLA4	ERG	FLT3	HIST2H3A	KIT	MTOR	PDCD1LG2	PTPRS	SHQ1	TFRC
ALOX12B	BTG1	CTNNA1	ERRF1	FLT4	HIST2H3C	KLF4	MUTYH	PDGFRA	PTPRT	SLIT2	TGFBR1
AMER1	BTK	CTNNB1	ESR1	FOXA1	HIST2H3D	KLHL6	MYB	PDGFRB	QKI	SLX4	TGFBR2
ANKRD11	CALR	CUL3	ETS1	FOXL2	HIST3H3	KMT2A	MYC	PDK1	RAB35	SMAD2	TMEM127
ANKRD26	CARD11	CUX1	ETV1	FOXO1	HNF1A	KRAS	MYCL	PDPK1	RAC1	SMAD3	TMPRSS2
APC	CASP8	CXCR4	ETV4	FOXP1	HNRNP	LAMP1	MYCN	PGR	RAD21	SMAD4	TNFAIP3
AR	CBFB	CYLD	ETV5	FRS2	HOXB13	LATS1	MYD88	PHF6	RAD50	SMARCA4	TNFRSF14
ARAF	CBL	DAXX	ETV6	FUBP1	HRAS	LATS2	MYOD1	PHOX2B	RAD51	SMARCB1	TOP1
ARFRP1	CCND1	DCUN1D1	EWSR1	FYN	HSD3B1	LMO1	NAB2	PIK3C2B	RAD51B	SMARCD1	TOP2A
ARID1A	CCND2	DDR2	EZH2	GABRA6	HSP90AA1	LRP1B	NBN	PIK3C2G	RAD51C	SMC1A	TP53
ARID1B	CCND3	DDX41	FAM46C	GATA1	ICOSLG	LYN	NCOA3	PIK3C3	RAD51D	SMC3	TP63
ARID2	CCNE1	DHX15	FANCA	GATA2	ID3	LZTR1	NCOR1	PIK3CA	RAD52	SMO	TRAF2
ARID5B	CD274	DICER1	FANCC	GATA3	IDH1	MAGI2	NEGR1	PIK3CB	RAD54L	SNCAIP	TRAF7
ASXL1	CD276	DIS3	FANCD2	GATA4	IDH2	MALT1	NF1	PIK3CD	RAF1	SOCS1	TSC1
ASXL2	CD74	DNAJB1	FANCE	GATA6	IFNGR1	MAP2K1	NF2	PIK3CG	RANBP2	SOX10	TSC2
ATM	CD79A	DNMT1	FANCF	GEN1	IGF1	MAP2K2	NFE2L2	PIK3R1	RARA	SOX17	TSHR
ATR	CD79B	DNMT3A	FANCG	GID4	IGF1R	MAP2K4	NFKBIA	PIK3R2	RASA1	SOX2	U2AF1
ATRX	CDC73	DNMT3B	FANCI	GLI1	IGF2	MAP3K1	NKX2-1	PIK3R3	RB1	SOX9	VEGFA
AURKA	CDH1	DOT1L	FANCL	GNA11	IKBKE	MAP3K13	NKX3-1	PIM1	RBM10	SPEN	VHL
AURKB	CDK12	E2F3	FAS	GNA13	IKZF1	MAP3K14	NOTCH1	PLCG2	RECQL4	SPOP	VTCN1
AXIN1	CDK4	EED	FAT1	GNAQ	IL10	MAP3K4	NOTCH2	PLK2	REL	SPTA1	WISP3
AXIN2	CDK6	EGFL7	FBXW7	GNAS	IL7R	MAPK1	NOTCH3	PMAIP1	RET	SRC	WT1
AXL	CDK8	EGFR	FGF1	GPS2	INHA	MAPK3	NOTCH4	PMS1	RHEB	SRSF2	XIAP
B2M	CDKN1A	EIF1AX	FGF2	GREM1	INHBA	MAX	NPM1	PMS2	RHOA	STAG1	XPO1
BAP1	CDKN1B	EIF4A2	FGF3	GRIN2A	INPP4A	MCL1	NRAS	PNRC1	RICTOR	STAG2	XRCC2
BARD1	CDKN2A	EIF4E	FGF4	GRM3	INPP4B	MDC1	NRG1	POLD1	RIT1	STAT3	YAP1
BBC3	CDKN2B	ELOC	FGF5	GSK3B	INSR	MDM2	NSD1	POLE	RNF43	STAT4	YES1
BCL10	CDKN2C	EML4	FGF6	H3F3A	IRF2	MDM4	NTRK1	PPARG	ROS1	STAT5A	ZBTB2
BCL2	CEBPA	EMSY	FGF7	H3F3B	IRF4	MED12	NTRK2	PPM1D	RPS6KA4	STAT5B	ZBTB7A
BCL2L1	CENPA	EP300	FGF8	H3F3C	IRS1	MEF2B	NTRK3	PPP2R1A	RPS6KB1	STK11	ZFHX3
BCL2L11	CHD2	EPCAM	FGF9	HGF	IRS2	MEN1	NUP93	PPP2R2A	RPS6KB2	STK40	ZNF217
BCL2L2	CHD4	EPHA3	FGF10	HIST1H1C	JAK1	MET	NUTM1	PPP6C	RPTOR	SUFU	ZNF703
BCL6	CHEK1	EPHA5	FGF14	HIST1H2BD	JAK2	MGA	PAK1	PRDM1	RUNX1	SUZ12	ZRSR2
BCOR	CHEK2	EPHA7	FGF19	HIST1H3A	JAK3	MITF	PAK3	PREX2	RUNX1T1	SYK	

Table 4: Additional content in TruSight Oncology Comprehensive (EU)

Fusions: 23 genes (from RNA)					
ALK	EGFR	ETV1	FGFR3	NTRK2	RET
AXL	EML4	ETV4	KIF5B	NTRK3	ROS1
BCL2	ERG	FGFR1	NRG1	PAX3	TMPRSS2
BRAF	ESR1	FGFR2	NTRK1	RAF1	
Splice variants: Two genes (from RNA)					
MET			EGFR		
Amplifications: Two genes (from DNA)					
ERBB2			MET		
Complex genomic signatures					
TMB			MSI		

What are the key attributes of the TruSight Oncology Comprehensive (EU) assay?

- Provides a kitted solution that can be implemented in-house by any lab
- Includes both DNA and RNA content and detects all variant classes, plus genomic signatures such as TMB and MSI; fusions are identified from RNA to maximize sensitivity for detection
- Enables CGP test results to be generated in only 4–5 days

How is TruSight Oncology Comprehensive (EU) different from current research use only (RUO) CGP assays, such as OncoPrint Comprehensive Assay Plus RNA*?

As a CE-marked IVD test, TruSight Oncology Comprehensive (EU) is compliant with European IVD Directive (IVDD) requirements and is on course to comply with stricter IVD Regulation (IVDR) legislation; the IVD label provides labs with the benefits of IVDR preparedness and reduced liability risk, and enables easier implementation with significantly reduced test validation efforts, as compared to RUO assays.

- Illumina has built a pipeline through multiple pharmaceutical partnerships based on TruSight Oncology Comprehensive (EU), aiming to expand CDx claims
- Built on proven Illumina technology for library preparation, sequencing, and bioinformatics, TruSight Oncology Comprehensive (EU) delivers highly reliable data with maximum quality and accuracy. Its hybrid-capture library preparation enables full detection and characterization of fusion events, which are not possible with amplicon-based techniques (Figure 1)




Figure 1: Hybrid-capture chemistry detects novel fusions missed by amplicon-based approaches

Amplicon-based approaches typically require confirmatory testing and do not characterize novel fusion partners. Hybrid-capture chemistry can identify both known and novel fusion partners.

* Thermo Fisher Scientific, Catalog no. A48578.

Workflow

What is the recommended sample input?

TruSight Oncology Comprehensive (EU) requires 40 ng RNA and/or 40 ng DNA extracted from FFPE tissue.

What are the workflow steps from sample preparation to final report?

The TruSight Oncology Comprehensive (EU) workflow includes four steps: sample acquisition and processing, DNA and RNA extraction, library preparation, and fully automated sequencing, analysis, and report generation (Figure 2).




Figure 2: TruSight Oncology Comprehensive (EU) sample-to-report workflow.

Additional details:

- Specimen tissue should be fixed using formalin fixative suitable for molecular analyses; a minimum of 20% tumor cell content is recommended and $\geq 30\%$ is optimal
- DNA and RNA extraction can be performed using commercially available extraction kits
- Library preparation takes approximately two days
- After loading the libraries on the NextSeq 550Dx instrument, the workflow is fully automated, including sequencing, base calling and QC, variant calling, interpretation, and generation of the final results report

How long is the turnaround time (TAT) from sample to report?

The TAT is 4–5 days from extracted DNA/RNA to the final results report.

What sequencing platform is needed?

TruSight Oncology Comprehensive (EU) is run on the NextSeq 550Dx instrument, an FDA-regulated and CE-marked, high-throughput sequencing platform.

What is the workflow for data analysis?

TruSight Oncology Comprehensive (EU) offers a streamlined, automated workflow from sequencing to final results report. Simply set up the sequencing run using Local Run Manager software. After the sequencing run is complete, secondary and tertiary analyses kickoff automatically and are performed on-instrument. The output is an easy-to-interpret results report.

What is the expected analysis time for a sample batch processed in a sequencing run?

The analysis time is 8–10 hours.

What is included in the final results report?

One of the key concerns when using a CGP panel is how to interpret the data and filter nonsignificant variants. TruSight Oncology Comprehensive (EU) software performs analysis and filtering as part of a fully automated workflow. The final report is easy to read and actionable.

Results Report

The Results Report is organized into two main sections ([Figure 3](#)):


1. Genomic findings with evidence of clinical significance:

Lists detected variants that have evidence of clinical significance (therapeutic, prognostic, or diagnostic) based on information in approved drug labels, guidelines, and clinical practice guidelines for the patient's tumor type.

2. Genomic findings with evidence of potential clinical significance:

Lists detected variants that:

- Have evidence of potential clinical significance (therapeutic, prognostic, or diagnostic) based on information in drug labels, guidelines, and clinical practice guidelines in another tumor type
- Match genomic and tumor type eligibility criteria for a clinical trial
- Have evidence of potential clinical significance in the primary literature for the patient's tumor type



Genomic findings with evidence of clinical significance

Genomic findings with evidence of potential clinical significance

Figure 3: TruSight Oncology Comprehensive (EU) example results report

The results report includes genomic findings with evidence of clinical significance and genomic findings with evidence of potential clinical significance.

Reimbursement

What reimbursement is available for CGP tests in Europe? How does that vary by country?

National and/or regional funding is available in most Western European countries. Laboratory services are reimbursed depending on the country in which they are located, the clinical setting in which they operate, and the services they provide (Figure 4 and Table 5). Illumina has established a dedicated Market Access team that works with Payers to expand the reimbursement of CGP across Europe. In addition, Illumina is working to open access to CGP in various major European emerging markets. For further questions about the appropriate coding or reimbursement for cancer testing, customers are encouraged to:

- Consult their respective national medical society for guidance regarding managing reimbursement claims and refer to any national or regional budget holders (payers) for guidance on coding and reimbursement (eg, national pathology societies)
- Contact their local Illumina representative, who can connect them with Illumina Market Access team. The team can provide guidance on where to find general reimbursement fee schedules, answer general reimbursement queries, or direct customers to key local institutions for further support†




Figure 4: CGP test coverage options across Europe

† As a manufacturer of IVD products, Illumina cannot advise customers on which specific reimbursement codes to use when billing for TruSight Oncology Comprehensive (EU).

Table 5: Laboratory fee schedules

Country	Relevant fee schedules
Austria	Not applicable
Belgium	National Institute for Health and Disability Insurance
Denmark	Not applicable
England	National Genomic Test Directory
Finland	Not applicable
France	L'Assurance Maladie billing and remuneration
Germany	KBV Uniform Evaluation Standard
Israel	National Health Basket
Italy ^a	Italy: Rates for Hospital and Outpatient Care
Netherlands ^b	Not applicable
Norway	Not applicable
Scotland	Not applicable
Spain ^a	Not applicable
Sweden	Not applicable
Switzerland	Federal Office of Public Health
<p>a. For Italy, the national fee schedule is informative only; for Italy and Spain, coverage is defined at the regional level. b. For the Netherlands, coverage is negotiated at the local level between hospitals and health insurers. Not applicable indicates that no fee schedule is available; diagnostics are funded through global budget allocation to laboratory customers.</p>	

Learn more →

[TruSight Oncology Comprehensive IVD Solutions](#)

1.800.809.4566 toll-free (US) | +1.858.202.4566 tel
techsupport@illumina.com | www.illumina.com

© 2025 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.