

Global Diversity Array-8 Manifest File Release Notes

V1.0 D1

April 19, 2021

FOR RESEARCH USE ONLY

© 2021 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CPro, CytoChip, DesignStudio, Epicentre, GAIx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Introduction

These release notes detail the key changes to the Global Diversity Array manifest files (.bpm &.csv formats) for v1.0 D1 since the v1.0 A1 initial release.

NEW FEATURES:

- None

IMPROVEMENTS:

- Optimized for CNV detection
- Validated indel and SNV calling against RefSeq
- Removed probes flagged for potential reduced performance by improved bioinformatics pipeline

ADDITIONAL SUPPORT FILES:

- List of 10,336 markers removed from manifest A1 and reason(s) for removal
- Updated annotation file
- Updated cluster file

Revision History

Revision	Date	Description of Change
v1.0 D1	19-APR-2021	Validated SNPs call with RefSeq and improved bioinformatics tool
v1.0 A1	14-OCT-2019	Initial release