

TruSight Tumor 15

Customer Release Notes

LRM Module v2.1.1

For MiSeq and MiniSeq systems

August 8, 2024

Introduction

These Release Notes detail the key features and known limitations for the TruSight Tumor 15 v2.1.1 Local Run Manager module.

NEW FEATURES:

N/A

FIXED ISSUES:

- Fixed issue where samples appeared duplicated when requeuing a run, with errors shown and no option to delete the duplicate samples.
- Fixed issue in the "Samples & Results" page of LRM where the option to view the results for each sample would not show.
- Fixed issue in the "Modules & Manifest" page of LRM where the download icon for the Report Definition file was only partially clickable.

Release History

Revision	Release Reference	Originator	Description of Change
00	CN 1109922	Joao Simao	Initial Release