



Ref. Certif. No.

US-TUVR-011077-A1

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME

CB TEST CERTIFICATE

Product

Next Generation Sequencer

Name and address of the applicant

Illumina, Inc.
5200 Illumina Way
San Diego CA 92122, USA

Name and address of the manufacturer

Illumina, Inc.
5200 Illumina Way
San Diego CA 92122, USA

Name and address of the factory

Illumina Singapore Pte. Ltd. North Tech Lobby 3 #02-13/18
29 Woodlands Industrial Park E1
Singapore 757716, Singapore

Ratings and principal characteristics

AC 100-240V, 50/60Hz; 80W;
Class I

Trademark (if any)

illumina

Customer's Testing Facility (CTF) Stage used

Model / Type Ref.

iSeq 100

Additional information (if necessary may also be reported on page 2)

Re-issuance of Certificate US-TUVR-11077 dated 3/16/2019 due to administrative modification (correction in test report). See associated test report for details.

A sample of the product was tested and found to be in conformity with

IEC 61010-1:2010
IEC 61010-2-010:2014
IEC 61010-2-081:2015
for national deviations see test report

As shown in the Test Report Ref. No. which forms part of this Certificate

31880498 009

This CB Test Certificate is issued by the National Certification Body



TÜVRheinland®

TUV Rheinland of North America, Inc.
12 Commerce Road
Newtown, CT 06470 - USA

Date: 24.10.2019

Signature:

Martin Glagla